



913 pH/DO Meter

914 pH/DO/Conductometer

Combined measuring instruments for dissolved oxygen, pH, and conductivity



# Handheld excellence for your most important quality parameters

The latest generation of our proven range of meters has become more powerful. Using the 913 pH/DO Meter or the 914 pH/DO/Conductometer, you can determine the pH value, concentration of dissolved oxygen, and conductivity fast, easily, and reliably. Count on trusted Swiss quality and enjoy the ease of use of our latest handheld meters.



## GLP-compliant documentation

Print out your measuring results directly in the lab, or save them as a CSV report.

- Flexible: the instrument supports thermoprinters, A4 printers, or impact printers
- Simple: CSV reports for data management

## Measure pH and dissolved oxygen in parallel

With the **913 pH/DO Meter**, you can measure the pH value and the concentration of dissolved oxygen simultaneously in the same sample.

- Reliable: the instrument informs you when the sensor cap needs to be exchanged
- Flexible: ideal for measurements in the field or in the laboratory



### Measure pH value, dissolved oxygen, and conductivity with the same instrument

With the **914 pH/DO/Conductometer**, you can measure all three parameters (pH, dissolved oxygen, and conductivity) on the same platform.


- Accurate: when measuring dissolved oxygen in solutions containing high salt concentrations, this can be compensated automatically
- Dependable: the instrument informs you when the sensor cap needs to be exchanged



### Optical measurement of dissolved oxygen

The new **O<sub>2</sub>-Lumitrode** sensor provides faster results for dissolved oxygen than conventional sensors.

- Fast: the O<sub>2</sub>-Lumitrode is the sensor for dissolved oxygen with the fastest response time on the market
- Maintenance-free: the O<sub>2</sub>-Lumitrode sensor is always ready to use
- Robust: the sensor membrane is protected from mechanical stress and interfering external light



### Reliable measuring results

The instrument informs you when sensor performance declines and the sensor cap needs to be replaced.

- Convenient: active monitoring of sensor performance
- Cost-optimized: there is no defined limit to the lifetime of the sensor. Instead, it depends on the intensity of usage
- Reliable: correct results over the entire lifetime of the sensor cap

## ORDERING INFORMATION

### 913 pH/DO Meter

2.913.0020      pH/DO Meter

2.913.0120      pH/DO Meter with Primatrode and O<sub>2</sub>-Lumitrode in a case with various accessories

2.913.0220      pH/DO Meter with stand

### 914 pH/DO/Conductometer

2.914.0030      pH/DO/Conductometer

2.914.0130      pH/DO/Conductometer with conductivity sensor and O<sub>2</sub>-Lumitrode in a carrying case with various accessories

2.914.0230      pH/DO/Conductometer with stand

### O<sub>2</sub>-Lumitrode

6.1116.000      O<sub>2</sub>-Lumitrode with calibration vessel and oxygen standard 0 %

6.2329.000      Oxygen standard 0 %

6.5623.000      Replacement set for O<sub>2</sub>-Lumitrode with calibration vessel and oxygen standard 0 %